



REFERENCE DATA

DWG. No.	
TITLE	Gas tank unit for Mini-CUBE series T-10-SUS outline

Mini-CUBE series



Flow sheet/Electric diagram



Close output, when below the set pressure.
(Non-voltage B contact output DC30V 0.5A or AC125V 0.5A or less)

Symbol	Name	Note	Q'ty
T	TANK	SUS 12.8L	1
PS	Pressure gauge with switch	050314	1
ST	Stop cock	100350	1

No.	NAME	SPECIFICATIONS	NOTES
SCALE	1 : 4	PROJECTION	
Rev.			
KOJIMA INSTRUMENTS INC.			
DESIGNED	DRAWN	CHECKED	APPROVED
DATE	3 / Mar. / 2015		