




No.	NAME	SPECIFICATIONS	NOTES
1	Flow meter	RK1450	
2	Flow controller	2203	
3	Ball valve		
4	Pressure gauge		
5	Pressure regulator	6600	
6	Gas outlet		
7	Gas inlet		

DWG. No.	KT-047032-E01
TITLE	Gas Blender GM-2B outline
CUSTOMER	



SCALE	1 : 4	PROJECTION	
Rev.			
KOJIMA INSTRUMENTS INC.			
DESIGNED	DRAWN	CHECKED	APPROVED
DATE	5 / April / 2004		