




| No. | NAME | SPECIFICATIONS | NOTES |
|-----|--------------------|----------------|-------|
| 1 | Flow meter | 2503F | |
| 2 | Pressure gauge | | |
| 3 | Pressure regulator | 6600 | |
| 4 | Gas outlet | | |
| 5 | Gas inlet | | |

| | |
|----------|----------------------------|
| DWG. No. | KT-127051-E01 |
| TITLE | Gas Blender PMG-1A outline |
| CUSTOMER | |



| | | | |
|-------------------------|-----------------|------------|----------|
| SCALE | 1 : 2.5 | PROJECTION | |
| Rev. | | | |
| KOJIMA INSTRUMENTS INC. | | | |
| DESIGNED | DRAWN | CHECKED | APPROVED |
| | | | |
| DATE | 23 / Aug / 2012 | | |