


No.	NAME	SPECIFICATIONS	NOTES
1	Touch panel		
2	Emergency stop switch		
3	Power switch		
4	Pressure regulator		
5	Gas inlet		
6	Gas inlet		
7	Gas inlet		
8	Gas inlet		
9	Gas outlet		
10	Input output terminal		
11	Circuit protector		
12	AC power inlet		
13	LAN terminal		
14	USB terminal		

DWG. No.	KT-147027-E01
TITLE	BRENDA BR-4C outline
CUSTOMER	



SCALE	1 : 4	PROJECTION	
Rev.			
KOJIMA INSTRUMENTS INC.			
DESIGNED	DRAWN	CHECKED	APPROVED
DATE	21 / April / 2015		