





No.	NAME	SPECIFICATIONS	NOTES
1	Touch panel		
2	Emergency stop switch		
3	Power switch		
4	Pressure regulator		
5	Gas inlet		
6	Gas inlet		
7	Gas outlet		
8	Input output terminal		
9	Circuit protector		
10	AC power inlet		
11	LAN terminal		
12	USB terminal		

DWG. No.	KT-157028-E01
TITLE	BRENDA BR-2C outline
CUSTOMER	



SCALE		1 : 4	PROJECTION	
Rev.				
KOJIMA INSTRUMENTS INC.				
DESIGNED		DRAWN	CHECKED	APPROVED
DATE		21 / April / 2015		