


No.	NAME	SPECIFICATIONS	NOTES
1	Touch panel		
2	Emergency stop switch		
3	Power switch		
4	Pressure regulator		
5	Gas inlet		
6	Gas inlet		
7	Gas inlet		
8	Gas outlet		
9	Input output terminal		
10	Circuit protector		
11	AC power inlet		
12	LAN terminal		
13	USB terminal		

DWG. No.	KT-157031-E01
TITLE	BRENDA BR-3C outline
CUSTOMER	



SCALE	1 : 4	PROJECTION	
Rev.			
KOJIMA INSTRUMENTS INC.			
DESIGNED	DRAWN	CHECKED	APPROVED
DATE	21 / April / 2015		